

CAT Hematoxylin

Counterstain
902-CATHE-112222

BIOCARE
M E D I C A L

Catalog Number: BRRCATHE-M, MM, GL

Description: 500 mL, 1000 mL, 1 gal

Intended Use:

For Research Use Only. Not for use in diagnostic procedures.

Summary & Explanation:

Biocare's CAT Hematoxylin is intended for use in the histologic demonstration of nuclear staining. Hematoxylin solutions mordanted with aluminum are routinely used in the histopathology laboratory for nuclear staining. This staining technique is used to make the critical distinction between a normal nucleus and an abnormal one. The traditional formulas for nuclear demonstration using an aluminum mordant are those of Delafield, Harris, Ehrich, Lillie and Mayer. All use aluminum alum (ammonium aluminum sulfate) as the mordant salt. Since the aluminum salts are not in themselves oxidizers, it is necessary to expose the hematoxylin solution to air or chemicals to affect the conversion of hematoxylin to hematein. These formulas are also used in special-staining procedures as counterstains for nuclear demonstration. The addition of acid to alum hematoxylin solutions is believed to increase the selectivity of the stain for the nuclei and to counteract the rapid oxidizing effects of chemical oxidizing agents. This latter function enables the solution to maintain some hematoxylin in equilibrium with the hematein to ensure a better stain. Glycerol tends to stabilize the system against over-oxidation and aids in preventing rapid evaporation.

Known Applications:

Any well-fixed tissue. Tissue can be frozen, paraffin-embedded, or celloidin-embedded.

Materials and Methods:

Reagents Provided:

Kit Catalog No.	Component Description	Quantity x Volume
BRRCATHE-M	CAT Hematoxylin	1 x 500 mL
BRRCATHE-MM	CAT Hematoxylin	1 x 1000 mL
BRRCATHE-GL	CAT Hematoxylin	1 x gallon

Supplied As:

Ready-to-use

Reconstitution, Dilution and Mixing:

CAT Hemoxlylin is provided ready-to-use. No reconstitution, dilution or mixing is required.

Materials and Reagents Required but Not Provided:

Microscope slides, positively charged
Desert Chamber* (Drying oven)
Positive and negative tissue controls
Xylene (Could be substituted with xylene substitute*)
Ethanol or reagent alcohol
Decloaking Chamber* (Pressure cooker)
Deionized or distilled water
Wash buffer*
Pretreatment reagents*
Enzyme digestion*
Avidin-Biotin Blocking Kit*(Labeled streptavidin kits only)
Peroxidase block*
Protein block*
Primary antibody*
Negative control reagents*
Detection kits*
Detection components*

Materials and Reagents Required but Not Provided Cont'd:

Chromogens*

Bleuing reagent*

Mounting medium*

* Refer to the Biocare Medical website located at <http://biocare.net> for information regarding catalog numbers and ordering.

Storage and Stability:

Store at room temperature away from light. The product is stable to the expiration date printed on the vial label when stored under these conditions. Do not use after expiration date. Storage under any condition other than those specified must be verified. The reagent(s) are ready-to-use and should not be diluted. The stability of user diluted reagent has not been established by Biocare.

Instructions for Use:

1. Deparaffinize and hydrate to water.
2. If sections are Zenker-fixed, remove the mercuric chloride crystals with iodine and clear with sodium thiosulfate (hypo). Rinse in running water.
3. Rinse in distilled water.
4. CAT Hematoxylin for 4-8 minutes. (*For IHC counterstain stain 30 seconds to 1 minute)
5. Wash in tap water until clear (generally 30 seconds to 1 minute).
6. Blue in Tacha's Bleuing Solution 30 seconds to 1 minute.
7. Wash in tap water 3-4 minutes.
8. Place in 95% alcohol for 30 seconds.
9. Counterstain with Edgar Degas Eosin 1-3 minutes. (*Alternately, counterstain with Rubens Eosin-Phloxine 30-40 seconds)
10. Dehydrate through graded alcohols.
11. Clear in clearing agents.
12. Mount with resinous mounting media.

Technical Notes:

Filter product with a 0.4-micron filter if particulate is seen.

Limitations:

This product is provided for Research Use Only (RUO) and is not for use in diagnostic procedures. Suitability for specific applications may vary and it is the responsibility of the end user to determine the appropriate application for its use.

Protocol Notes:

Biocare Medical uses only stains that are certified by the Biological Stain Commission. CAT Hematoxylin is specially formulated to eliminate the necessity for differentiation of the section.

NOTE: The biggest objection to Mayer's hematoxylin as used in the past, has been that stained slides often fade after 1-3 years. This problem can be eliminated, however, when slides are washed after hematoxylin in running water for a minimum of 3-4 minutes.

NOTE: Biocare's research studies show pH of tap and DI water for rinsing slides can influence staining. pH can change from acidic to basic seasonally in different regions of the country. It has been further reported that there have been changes in tap water from the morning to the afternoon in some states. Long rinsing times after bleuing may differentiate out Hematoxylin. We have modified our rinse time to 3-4 minutes, instead of the previous 7-10 minutes.

Precautions:

1. Refer to reagent Safety Data Sheet for precautions.
2. Contains glacial acetic acid (use in a well-ventilated areas).
3. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting



60 Berry Drive
Pacheco, CA 94553

USA

Rev. 062117

Tel: 800-799-9499 | www.biocare.net | Fax: 925-603-8080

CAT Hematoxylin

Counterstain
902-CATHE-112222


BIOCARE
M E D I C A L

Precautions Cont'd:

the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.¹

4. Microbial contamination of reagents may result in an increase in nonspecific staining.
5. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.
6. Do not use reagent after the expiration date printed on the vial.

This CAT Hematoxylin contains components classified as indicated in the table below in accordance with the Regulation (EC) No. 1272/2008

Hazard	Code	Hazard Statement
	H315 H319	Causes skin irritation. Causes serious eye irritation.

Technical Support:

Contact Biocare's Technical Support at 1-800-542-2002 for questions regarding this product.

References:

1. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition (M29-A4) Wayne, PA 2014.