p16 + Ki-67

Prediluted Multiplex Antibody Reagent
901-3246DS-071223

English

BIOCARE
M E D I C A L

Available Product Formats				
Format	Catalog Number	Description	Dilution	Diluent
ONCORE Pro	OPAI 3246DS T60	60 tests	Ready-to-use	N/A

Intended Use:

For In Vitro Diagnostic Use

p16 + Ki-67 is a cocktail of mouse monoclonal and rabbit monoclonal antibodies that is intended for laboratory use in the qualitative identification of p16 and Ki-67 proteins by immunohistochemistry (IHC) in formalin-fixed paraffin-embedded (FFPE) human tissues. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation:

p16 INK4a is a tumor suppressor protein involved in the pathogenesis of a variety of malignancies. p16 INK4a is a 16.5 kDa protein expressed in the nucleoplasm of proliferating cells, functioning as an inhibitor of CDK4. Recent analyses of the p16 INK4a gene revealed homozygous deletions, nonsense, missense, or frameshift mutations in several human cancers.¹

The Ki-67 nuclear antigen is associated with cell proliferation. It is found throughout the cell cycle in the G1, S, G2, and M phases; but not the (G0) phase. It is commonly used to grade proliferation rates of tumors.

An optimized antibody cocktail for p16 INK4a and Ki-67 may aid in the identification of cells co-expressing markers for both tumor suppression and cell proliferation, an indicator of cell-cycle deregulation.²

Principle of Procedure:

Antigen detection in tissues and cells is a multi-step immunohistochemical process. The initial step binds the primary antibody to its specific epitope. After labeling the antigen with a primary antibody, a one-step or two-step detection procedure can be applied. A one-step procedure will feature an enzyme labeled polymer that binds the primary antibody. A two-step procedure will feature a linker antibody added to bind to the primary antibody. An enzyme-labeled polymer is then added to bind the linker antibody. These detections of the bound antibodies are evidenced by a colorimetric reaction.

Reagent Provided:

p16 + Ki-67 is provided as a prediluted antibody cocktail of anti-p16 and anti-Ki-67 antibodies in buffer with carrier protein and preservative.

Antibody	anti-p16	anti-Ki-67
Clone	BC42	SP6
Source	Mouse monoclonal	Rabbit monoclonal
Isotype	IgG1/kappa	IgG
Epitope/ Antigen	p16 INK4a	Ki-67
Cellular Localization	Nuclear and cytoplasmic	Nuclear
Staining	Brown (DAB)	Red

Species Reactivity: Human; others not tested

Protein Concentration: Call for lot specific Ig concentration.

Positive Tissue Control: Cervical Cancer

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Storage and Stability:

Store at 2°C to 8°C. The product is stable to the expiration date printed on the label, when stored under these conditions. Do not use after expiration date. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

Protocol Recommendations (ONCORE™ Pro Automated Slide Staining System):

OPAI3246DS is intended for use with the ONCORE Pro. Refer to the User Manual for specific instructions for use. Protocol parameters in the Protocol Editor should be programmed as follows:

Protocol Name: p16 + Ki-67

Protocol Template (Description): Multiplex 2 Template 1

Dewaxing (DS Option): DS2-50

Antigen Retrieval (AR Option): AR1, high pH; 103°C

Block Option: Buffer

Reagent Name, Time, Temp.: p16 + Ki-67, 59 min., 25°C

ONCORE Pro Template Options	
Template 1 (Default)	Default Templates (No Hematoxylin, AP protocols utilize Fast Red)
Template 2	AP & Multiplex Templates with Warp Red
Template 3	Default HRP Templates + Hematoxylin
	AP & Multiplex Templates with Warp Red + Hematoxylin

Limitations:

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to fixation, heat-retrieval method, incubation times, tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of Biocare products. Ultimately, it is the responsibility of the investigator to determine optimal conditions.

Quality Control:

Refer to CLSI Quality Standards for Design and Implementation of Immunohistochemistry Assays; Approved Guideline-Second edition (I/LA28-A2) CLSI Wayne, PA USA (www.clsi.org). 2011

Precautions:

1. This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard communication and EC Directive 91/155/EC. Sodium azide (NaN₃) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for Disease Control, 1976, National Institute of Occupational Safety and Health, 1976)³
2. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come into contact with sensitive areas, wash with copious amounts of water.⁴
3. Microbial contamination of reagents may result in an increase in nonspecific staining.
4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.
5. Do not use reagent after the expiration date printed on the vial.
6. The SDS is available upon request and is located at <http://biocare.net>.

Troubleshooting:

Follow the antibody specific protocol recommendations according to data sheet provided. If atypical results occur, contact Biocare's Technical Support at 1-800-542-2002.

p16 + Ki-67

Prediluted Multiplex Antibody Reagent

901-3246DS-071223

English

BIOCARE
M E D I C A L

References:

1. LaPak KM, Burd CE. The molecular balancing act of p16(INK4a) in cancer and aging. *Mol Cancer Res.* 2014 Feb;12(2):167-83.
2. Ikenberg H, et al. Screening for cervical cancer precursors with p16/Ki-67 dual-stained cytology: results of the PALMS study. *J Natl Cancer Inst.* 2013;105:1550-7.
3. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
4. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

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