

CDX2 + CK7

Prediluted Multiplex Cocktail (4-Step)

Control Number: 901-367DS-090817

Catalog Number: PM 367 DS AA, H, L
Description: 6.0, 25, 100ml, prediluted
Dilution: Ready-to-use
Diluent: N/A

Intended Use:
 For In Vitro Diagnostic Use

Summary and Explanation:

CDX2 is a homeobox gene that encodes an intestine-specific transcription factor. It is expressed in the nuclei of epithelial cells throughout the intestine, from duodenum to rectum. The CDX2 protein is expressed in primary and metastatic colorectal carcinomas and has also been demonstrated in the intestinal metaplasia of the stomach and intestinal-type gastric cancer, while it is not expressed in the normal gastric mucosa. Studies have shown that CDX2 is a superior marker compared to CK20 and can be substituted in a panel of antibodies.

Cytokeratin 7 is a basic cytokeatin and is expressed in epithelial cells of ovary, lung and breast, but not in the colon or gastrointestinal tract. It is often used in conjunction with Cytokeratin 20 in distinguishing pulmonary ovarian and breast carcinomas (CK7 +) from colon carcinomas (CK7-).

This Multiplex cocktail of CDX2 and CK7 can be used to distinguish colon cancers from breast, lung and ovarian cancers. CDX-2 will stain the nuclei brown and CK7 will stain target antigens red.

Principle of Multiplex Staining:

A Multiplex IHC stain can be accomplished in four major steps. The initial step consists of an antibody cocktail with at least one mouse and one rabbit antibody. This cocktail is applied to the tissue and will bind with two or more target antigens. A multiplex detection cocktail of horseradish peroxidase (HRP) and alkaline phosphatase (AP) conjugated secondary antibodies is applied. The third step consists of the addition of DAB-Substrate that binds to the HRP and produces a brown chromogenic reaction product. The fourth step consists of a Fast Red-Substrate that binds to the AP and produces a red chromogenic reaction product.

Source: Mouse Monoclonal and Rabbit Monoclonal

Species Reactivity: Human; others not tested

Clone: CDX2-88 + BC1

Isotype: IgG1 and Rabbit IgG

Epitope/Antigen: CDX2 and CK7

Cellular Localization:

CDX2: (nuclear); brown

CK7: (cytoplasmic); red

Positive Control: Colon, breast, ovary and lung cancers

Normal Tissue: Colon, breast

Abnormal Tissue: Colon, breast and lung cancer

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Supplied As: Buffer with protein carrier and preservative.

Storage and Stability:

Store at 2°C to 8°C. Do not use after expiration date printed on vial. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

Protocol Recommendations

Peroxide Block:

Block for 5 minutes with BIOCARE's PEROXIDAZED 1.

Pretreatment Solution (recommended): Diva

Pretreatment Protocol:

Heat Retrieval Method:

Retrieve sections under pressure using BIOCARE's Decloaking Chamber, followed by a wash in distilled water. Alternatively, steam tissue sections for 45-60 minutes. Allow solution to cool for 20 minutes then wash in distilled water.

Protein Block:

Optional: Incubate for 10-15 minutes at RT with BIOCARE's Background Sniper.

Primary Antibody:

Incubate for 30 minutes at RT.

Double Stain Detection:

Incubate for 30 minutes at RT using BIOCARE's MACH 2 Double Stain 2.

Chromogen (1): Incubate for 5 minutes at RT when using BIOCARE's Betazoid DAB.

Chromogen (2):

Incubate for 10-20 minutes at RT with BIOCARE's Vulcan Fast Red. Rinse in deionized water.

Technical Notes:

This antibody has been standardized with BIOCARE's MACH 2 Double Stain 2. It can also be used on an automated staining system. Use TBS buffer for washing steps.

Limitations:

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to: fixation, heat-retrieval method, incubation times, tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of BIOCARE products. Ultimately, it is the responsibility of the investigator to determine optimal conditions. These products are tools that can be used for interpretation of morphological findings in conjunction with other diagnostic tests and pertinent clinical data by a qualified pathologist.

Quality Control:

Refer to NCCLS Quality Assurance for Immunocytochemistry approved guidelines, December 1999 MM4-A Vol.19 No.26 for more information about Tissue Controls.

Precautions:

This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard communication and EC Directive 91/155/EC.

Sodium azide (NaN₃) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for disease control, 1976, National Institute of Occupational Safety and Health, 1976)

Specimens, before and after fixation and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Microbial contamination of reagents may result in an increase in nonspecific staining. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change. The MSDS is available upon request.

Troubleshooting:

Follow the antibody specific protocol recommendations according to data sheet provided. If atypical results occur, contact BIOCARE's Technical Support at 1-800-542-2002.

CDX2 + CK7

Prediluted Multiplex Cocktail (4-Step)

Control Number: 901-367DS-090817

References:

1. Kim MJ. The usefulness of CDX-2 for differentiating primary and metastatic ovarian carcinoma: an immunohistochemical study using a tissue microarray. *J Korean Med Sci.* 2005 Aug;20(4):643-8.
2. Kennedy MT, et al. Expression pattern of CK7, CK20, CDX-2, and villin in intestinal-type sinonasal adenocarcinoma. *J Clin Pathol.* 2004 Sep;57(9):932-7.
3. Werling RW, Yaziji H, Bacchi CE, Gown AM. CDX2, a highly sensitive and specific marker of adenocarcinomas of intestinal origin: an immunohistochemical survey of 476 primary and metastatic carcinomas. *Am J Surg Pathol.* 2003 Mar;27(3):303-10.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
5. National Committee for Clinical Laboratory Standards(NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, PA 1991;7(9). Order code M29-P.