## Bcl-6 [LN22]

Prediluted Monoclonal Antibody  
901-410-081324

English

**BIOCARE**  
M E D I C A L

### Available Product Formats

Format	Catalog Number	Description	Dilution	Diluent
Q Series— For Leica BOND-III	ALI 410 G7	7.0 mL	Ready-to-use	N/A

#### Intended Use:

For In Vitro Diagnostic Use

Bcl-6 [LN22] is a mouse monoclonal antibody that is intended for laboratory use in the qualitative identification of Bcl-6 protein by immunohistochemistry (IHC) in formalin-fixed paraffin-embedded (FFPE) human tissues. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Summary and Explanation:

The Bcl-6 gene, originally cloned from a tumor with 3q27 translocation, is commonly expressed in large cell lymphomas. In humans, Bcl-6 encodes for a Kruppel-type zinc finger protein that is believed to be important in germinal center formation. Bcl-6 protein is expressed mainly by follicle center cells and in a few interfollicular T-lymphocytes. Bcl-6 has also been detected in Burkitt's lymphoma and in nodular lymphocyte predominant Hodgkin's disease. Bcl-6 is not expressed in hairy cell leukemia, mantle cell or marginal-zone derived lymphomas. This antibody demonstrates strong nuclear reactivity.

#### Principle of Procedure:

Antigen detection in tissues and cells is a multi-step immunohistochemical process. The initial step binds the primary antibody to its specific epitope. After labeling the antigen with a primary antibody, a one-step or two-step detection procedure can be applied. A one-step procedure will feature an enzyme labeled polymer that binds the primary antibody. A two-step procedure will feature a linker antibody added to bind to the primary antibody. An enzyme-labeled polymer is then added to bind the linker antibody. These detections of the bound antibodies are evidenced by a colorimetric reaction.

**Source:** Mouse monoclonal

**Species Reactivity:** Human

**Clone:** LN22

**Isotype:** IgG2b

**Protein Concentration:** Call for lot specific Ig concentration.

**Epitope/Antigen:** AA 1-350 N-terminus human Bcl-6 oncoprotein

**Cellular Localization:** Nuclear

**Positive Tissue Control:** Tonsil or Follicular lymphoma

**Known Applications:**

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

**Supplied As:** Buffer with protein carrier and preservative

**Storage and Stability:**

Store at 2°C to 8°C. The product is stable to the expiration date printed on the label, when stored under these conditions. Do not use after expiration date. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

#### Protocol Recommendations (Q Series – For Leica BOND-III):

ALI410 is intended for use with the Leica BOND-III. Refer to the User Manual for specific instructions for use. Recommended protocol parameters are as follows:

**Protocol Name:** IHC Protocol F

**Detection:** Bond polymer Refine

**HIER:** 40 min with ER2

**Peroxide Block:** 5 min

**Marker (Primary Antibody):** 15 min

**Post Primary:** 8 min

**Polymer:** 8 min

**Mixed DAB Refine:** 10 min

**Hematoxylin:** 5 min

#### Limitations:

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to fixation, heat-retrieval method, incubation times, tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of Biocare products. Ultimately, it is the responsibility of the investigator to determine optimal conditions.

#### Quality Control:

Refer to CLSI Quality Standards for Design and Implementation of Immunohistochemistry Assays; Approved Guideline-Second edition (I/LA28-A2) CLSI Wayne, PA USA ([www.clsi.org](http://www.clsi.org)). 2011

#### Precautions:

1. This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard communication and EC Directive 91/155/EC. Sodium azide (NaN<sub>3</sub>) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for Disease Control, 1976, National Institute of Occupational Safety and Health, 1976)<sup>4</sup>

2. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come into contact with sensitive areas, wash with copious amounts of water.<sup>5</sup>

3. Microbial contamination of reagents may result in an increase in nonspecific staining.

4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.

5. Do not use reagent after the expiration date printed on the vial.

6. The SDS is available upon request and is located at <http://biocare.net>.

#### Troubleshooting:

Follow the antibody specific protocol recommendations according to data sheet provided. If atypical results occur, contact Biocare's Technical Support at 1-800-542-2002.

#### References:

1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. *Leuk Lymphoma*. 2001 May; 41(5-6):585-92.

2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. *Am J Surg Pathol*. 2001 Mar; 25(3):307-15.

3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. *Am J Surg Pathol*. 2000 May; 24(5):694-702.

4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."

5. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

Q Series antibodies are developed solely by Biocare Medical LLC and do not imply approval or endorsement of Biocare antibodies by Leica Biosystems. Biocare and Leica Biosystems are not affiliated, associated or related in any way. Leica, Leica Biosystems, BOND-MAX and BOND-III are trademarks of Leica Biosystems.

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1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. Leuk Lymphoma. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. Am J Surg Pathol. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. Am J Surg Pathol. 2000 May;24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."

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1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. *Leuk Lymphoma*. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. *Am J Surg Pathol*. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. *Am J Surg Pathol*. 2000 May; 24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."

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2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. Am J Surg Pathol. 2001 Mar; 25(3):307-15.

3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. Am J Surg Pathol. 2000 May; 24(5):694-702.

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#### Referencer:

1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. Leuk Lymphoma. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. Am J Surg Pathol. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. Am J Surg Pathol. 2000 May; 24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
5. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

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#### Références :

1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. *Leuk Lymphoma*. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. *Am J Surg Pathol*. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. *Am J Surg Pathol*. 2000 May; 24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
5. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

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#### Referenzen:

1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. Leuk Lymphoma. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. Am J Surg Pathol. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. Am J Surg Pathol. 2000 May; 24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."

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- Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. *Leuk Lymphoma*. 2001 May; 41(5-6):585-92.
- Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. *Am J Surg Pathol*. 2001 Mar; 25(3):307-15.
- Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. *Am J Surg Pathol*. 2000 May; 24(5):694-702.
- Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
- Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

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#### Referencias:

1. Dunphy CH, *et al.* Paraffin immunoreactivity of CD10, CDw75, and Bcl-6 in follicle center cell lymphoma. *Leuk Lymphoma*. 2001 May; 41(5-6):585-92.
2. Fernández-Vázquez A, *et al.* Primary cutaneous large B-cell lymphoma: the relation between morphology, clinical presentation, immunohistochemical markers, and survival. *Am J Surg Pathol*. 2001 Mar; 25(3):307-15.
3. Yang B, *et al.* Clinicopathologic reassessment of primary cutaneous B-cell lymphomas with immunophenotypic and molecular genetic characterization. *Am J Surg Pathol*. 2000 May;24(5):694-702.
4. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."
5. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

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