

Val Wash Buffer, 20X

Wash Buffer
901-VLT8013-112718



BIOCARE
M E D I C A L

Catalog Number: VLT 8013 MX
Description: 500 mL, concentrated

Intended Use:

For In Vitro Diagnostic Use
Val Wash Buffer, 20X (sodium azide and thimerosal free) is intended for use in automated immunohistochemistry applications on the VALENT® Automated Slide Staining Platform. The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary & Explanation:

Val Wash Buffer, 20X (sodium azide and thimerosal free) is a high quality TBS buffer with superior pH stability. The buffer contains a surfactant that allows reagents to spread uniformly across tissue sections on slides. It is designed for use on the VALENT Automated Slide Staining Platform.

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Reagents Provided:

Val Wash Buffer, 20X (VLT8013MX) 1 x 500 mL

Reconstitution, Dilution and Mixing:

Val Wash Buffer, 20X is provided as a concentrate. See Instructions for Use for proper dilution.

Materials and Reagents Required But Not Provided:

Reagents and materials, such as primary antibodies, detection kits, chromogens and ancillary reagents are not provided.

Refer to the Biocare Medical website located at <http://biocare.net> for information regarding catalog numbers and ordering information.

Refer to the VALENT Automated Slide Staining Platform User Manual for a complete list of VALENT specific materials and reagents required.

Storage and Stability:

Store at room temperature. The product is stable to the expiration date printed on the label, when stored under these conditions. Do not use after expiration date. Diluted reagents should be used promptly; any remaining reagent should be stored at room temperature.

Instructions for Use:

1. Mix 1 part concentrate to 19 parts deionized water (1:20 dilution), or dilute contents of the Val Wash Buffer, 20X bottle (500 mL) to 9.5 liters of deionized water.
2. Check pH. If necessary, adjust to 7.7 ± 0.1 at 25°C.

Val Wash Buffer, 20X is provided as a concentrate for the VALENT Automated Slide Staining Platform. Place the ready-to-use diluted Val Wash Buffer in the appropriately labeled Valent bulk carboy. The VALENT Automated Slide Staining Platform will apply reagent as required in the selected protocol.

Refer to the appropriate antibody data sheet for the recommended staining protocol. Refer to the VALENT Automated Slide Staining Platform User Manual for detailed instructions on instrument operation and additional protocol options.

Technical Notes:

The diluted Val Wash Buffer is stable for 1 week. The carboy used for the Val Wash Buffer should be cleaned and sanitized once per month. Clean with 1 liter of 70% alcohol. Wash well with deionized water. Failure to sanitize the carboy may result in bacterial and mold growth in the carboy.

Limitations:

These reagents have been optimized for use with VALENT antibodies and ancillary reagents. The protocols for a specific application can vary. These include, but are not limited to fixation, heat-retrieval method, incubation times, and tissue section thickness. Third party primary antibodies may be used on the VALENT Automated Slide Staining Platform; however, appropriate antibody concentration may depend upon multiple factors and must be empirically determined by the user. Ultimately, it is the responsibility of the investigator to determine optimal conditions.

Quality Control:

Refer to CLSI Quality Standards for Design and Implementation of Immunohistochemistry Assays; Approved Guideline-Second edition (I/LA28-A2) CLSI Wayne, PA USA (www.clsi.org). 2011

Precautions:

1. Refer to reagent Safety Data Sheet for precautions.
2. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. (1)
3. Microbial contamination of reagents may result in an increase in nonspecific staining.
4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.
5. Do not use reagent after the expiration date printed on the vial.




IRRITANT

Troubleshooting:

Follow the antibody specific protocol recommendations according to data sheet provided. If atypical results occur, contact Biocare's Technical Support at 1-800-542-2002.

References:


1. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

 Biocare Medical
60 Berry Drive
Pacheco, CA 94553
USA



Rev: 073018

Tel: 800-799-9499 | www.biocare.net | Fax: 925-603-8080

 EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands