

# TTF-1 [SPT24] Staining Specificity in Lung Adenocarcinoma vs. Lung Squamous Cell Carcinoma is Markedly Improved with Titer Optimization and Cut-Off Values

David Tacha, PhD and Ding Zhou, BS; Biocare Medical, Concord, CA

---

## Introduction

Lung cancer is the leading cause of cancer death for both men and women. More people die of lung cancer than of colon, breast and prostate cancers combined.<sup>1</sup> Classification of lung carcinomas into histological types is typically performed by visual examination using Hematoxylin and Eosin (H&E) and immunohistochemistry (IHC). However, accurate classification can be difficult with poorly differentiated or undifferentiated lung carcinoma. Historically, antibodies TTF-1 and p63 have often been used to differentiate primary lung adenocarcinoma (LADC) from squamous cell carcinoma (SqCC) of the lung.<sup>2-4</sup>

Commercially available Thyroid Transcription Factor-1 (TTF-1) monoclonal antibodies 8G7G3/1 and SPT24 have been shown to have different sensitivities in LADC and in lung SqCC.<sup>5-7</sup> A study by Masai *et al.* demonstrated that SPT24 was more sensitive than 8G7G3/1 in LADC (72.4% and 65.4%, respectively).<sup>5</sup> However, the study demonstrated that SPT24 stained a higher percentage of lung SqCC (16.8% vs. 1%). Increased staining of SPT24 in lung SqCC has also been shown to be heavily influenced by different detection systems.<sup>8</sup> These findings with SPT24 led us to investigate antibody dilution factor and staining cut-off values along with additional screening antibodies to provide better specificity. Desmoglein 3 (DSG3) and p40 have recently been shown to be highly specific for lung SqCC while Napsin A has demonstrated absolute specificities in LADC.<sup>4,9,10</sup> Double stain cocktails including DSG3 + p40 as well as TTF-1 + Napsin A have shown high specificity and sensitivity in lung SqCC.<sup>4,10,11</sup> The co-expression of TTF-1 and Napsin A has also been shown to be highly specific for LADC, and its co-expression is regarded as pulmonary specific.<sup>11</sup>

In this study, we will screen cases of lung SqCC and LADC with a double stain cocktail of DSG3/p40 + Napsin A. The SPT24 titer will be optimized for LADC and cut-off values will be examined to improve overall staining specificity without compromising staining sensitivity.

## Materials and Methods

Formalin-fixed paraffin-embedded lung SqCC (n=137) and LADC (n=60) tissue microarrays were acquired. DSG3/p40 + Napsin A double stain cocktail (Biocare Medical, Concord, CA), TTF-1 [SPT24] and TTF-1 [8G7G3/1] antibody titers were optimized for IHC. SPT24 was optimally titered at 1:1200. The double stain antibody cocktail was detected by using a polymer double stain detection kit and visualization was achieved with DAB and Fast Red chromogens.


## Results

Results are summarized in Tables 1 and 2. TTF-1 [SPT24] and TTF-1 [8G7G3/1] both stained 1.5% (2/137) of lung SqCC (Table 1). Napsin A was negative in all lung SqCC while DSG3/p40 stained 91.2% of lung SqCC (Table 1). In LADC, TTF-1 [SPT24] stained 88.3% of cases compared to Napsin A (73.3%) and TTF-1 [8G7G3/1] (63.3%) (Table 2). Lung SqCC cases that demonstrated expression of SPT24 below the cutoff value (<10%) were positive for DSG3 and/or p40.


## Scoring and Interpretation Method

TTF-1 cases were considered positive if 10% or more of tumor cells were stained with a staining intensity of > 1+. Cases with <10% staining and no focal areas of positive staining were scored as negative.


Figure 1



Lung adenocarcinoma (LADC) stained with SPT24




LADC stained with 8G7G3/1




LADC stained with Napsin A


Figure 2




Liver stained with 8G7G3/1




Liver stained with SPT24 1:200




Liver stained with SPT24 1:1200



Lung SqCC stained with 8G7G3/1



Lung SqCC stained with SPT24 1:1200



Lung SqCC stained with SPT24 1:200

Figure 3




Table 1

DSG3/p40	Napsin A	TTF-1 [SPT24]	TTF-1 [8G7G3/1]
125/137 (91.2%)	0/137 (0%)	2/137 (1.5%)	2/137 (1.5%)

Table 2

DSG3/p40	Napsin A	TTF-1 [SPT24]	TTF-1 [8G7G3/1]
1/60 (1.7%)*	44/60 (73.3%)	53/60 (88.3%)	38/60 (63.3%)

\*One case of lung adenocarcinoma stained with p40 but also stained with Napsin A.

## Discussion

The well published use of 8G7G3/1 has been shown to stain lung adenocarcinoma and small cell lung carcinoma but, in most cases, not lung squamous cell carcinoma.<sup>2,3</sup> In this study, SPT24 and Napsin A were shown to be superior at detecting LADC when compared to 8G7G3/1 (Table 2, Figures 1A-C). When applying a wide range of SPT24 antibody titers (1:200-1:1200) on liver and lung cancer, we observed staining in normal liver (cytoplasmic) and in certain cases of lung SqCC when a titer of 1:200 was used. Normal liver and the expression of SPT24 originally observed in lung SqCC were both negative when SPT24 was systematically diluted out to 1:1200 (Figures 2A-F). Only residual normal lung was stained with both TTF-1 antibodies (Arrows, Figure 2D, 2E). Strong and diffuse staining in LADC was achieved with a 1:1200 dilution of SPT24 and staining intensity was comparable to 8G7G3/1 (Figures 3A, 3B). However, at 1:200, SPT24 staining intensity was stronger than the optimized titers of 8G7G3/1 and SPT24 (Figure 3C).

Based on the interpretation method in a study by Mukhopadhyay and Katzenstein, a  $\geq 10\%$  cut-off value for positive staining was used to increase specificity and decrease potential false positives.<sup>12</sup> Correlation of confirmed cases of LADC and lung SqCC was achieved by IHC screening with the double stain cocktail of DSG3/p40 + Napsin A (Tables 1 and 2). Two poorly differentiated cases that had been previously diagnosed as lung SqCC expressed both SPT24 and 8G7G3/1 antibodies but were negative for DSG3, p40 and Napsin A. As a result, we could not confirm the previous diagnosis of lung SqCC by IHC or by morphological examination.

## Conclusion

TTF-1 [SPT24] was more sensitive than TTF-1 [8G7G3/1] in LADC. TTF-1 [SPT24] specificity in lung SqCC is significantly improved with titer optimization and cut-off values and has been shown to be equivalent to TTF-1 [8G7G3/1]. When properly titered, TTF-1 [SPT24] does not stain normal liver whereas TTF-1 [8G7G3/1] does stain normal liver. The co-expression of TTF-1 and Napsin A provides further specificity for lung adenocarcinoma or when determining tumors of unknown origin. DSG3 and p40 are highly specific for lung SqCC and therefore, may be used as a pre-screener to rule out LADC. Caution or the use of other confirmatory markers should be considered when poorly differentiated lung cancers are expressing TTF-1 but are negative for DSG3, p40 and Napsin A.

## References

1. American Cancer Society: <http://www.cancer.org/cancer/lungcancer/index>.
2. Wu M, *et al.* p63 and TTF-1 immunostaining. A useful marker panel for distinguishing small cell carcinoma of lung from poorly differentiated squamous cell carcinoma of lung. *Am J Clin Pathol.* 2003 May;119(5):696-702.
3. Di Loreto C, *et al.* Immunocytochemical expression of tissue specific transcription factor-1 in lung carcinoma. *J Clin Pathol.* 1997 Jan;50(1):30-2.
4. Tacha D, *et al.* A 6-antibody panel for the classification of lung adenocarcinoma versus squamous cell carcinoma. *Appl Immunohistochem Mol Morphol.* 2012 May;20(3):201-7.
5. Masai K, *et al.* Expression of squamous cell carcinoma markers and adenocarcinoma markers in primary pulmonary neuroendocrine carcinomas. *Appl Immunohistochem Mol Morphol.* 2013 Jul;21(4):292-7.
6. Matoso A, *et al.* Comparison of thyroid transcription factor-1 expression by 2 monoclonal antibodies in pulmonary and nonpulmonary primary tumors. *Appl Immunohistochem Mol Morphol.* 2010 Mar;18(2):142-9.
7. Ordóñez NG. Value of thyroid transcription factor-1 immunostaining in tumor diagnosis: a review and update. *Appl Immunohistochem Mol Morphol.* 2012 Oct;20(5):429-44.
8. Kashima K, *et al.* Significant expression of thyroid transcription factor-1 in pulmonary squamous cell carcinoma detected by SPT24 monoclonal antibody and CSA-II system. *Appl Immunohistochem Mol Morphol.* 2014;22(2):119-24.
9. Tacha D, *et al.* An immunohistochemical analysis of a newly developed mouse monoclonal p40 (BC28) antibody in lung, bladder, skin, breast, prostate, and head and neck cancers. *Arch Pathol Lab Med.* 2014 Oct;138(10):1358-64.
10. Brown AF, *et al.* Tissue-preserving antibody cocktails to differentiate primary squamous cell carcinoma, adenocarcinoma, and small cell carcinoma of lung. *Arch Pathol Lab Med.* 2013 Sep;137(9):1274-81.
11. Ye J, *et al.* Combination of napsin A and TTF-1 immunohistochemistry helps in differentiating primary lung adenocarcinoma from metastatic carcinoma in the lung. *Appl Immunohistochem Mol Morphol.* 2011 Jul;19(4):313-7.
12. Mukhopadhyay S, Katzenstein AL. Subclassification of non-small cell lung carcinomas lacking morphologic differentiation on biopsy specimens: Utility of an immunohistochemical panel containing TTF-1, napsin A, p63, and CK5/6. *Am J Surg Pathol.* 2011 Jan;35(1):15-25.

**BIOCARE**  
M E D I C A L

800.799.9499  
4040 Pike Lane  
Concord CA 94520

[www.biocare.net](http://www.biocare.net)