

Anti-TIGIT Rabbit Monoclonal Antibody
[BLR047F] - A new member of the
immune checkpoint receptors family

Anti-TIGIT Rabbit Monoclonal Antibody [BLR047F]

A New Member of the Immune Checkpoint Receptors Family

Antibodies targeting immune checkpoint receptors, like cytotoxic T lymphocyte antigen-4 (CTLA-4) and cell death protein-1 (PD-1), have recently been utilized in novel immunotherapies as they have demonstrated remarkable clinical efficiency in different tumor types, such as metastatic melanoma, lung cancer, renal, and bladder carcinoma.¹ Because of its success, blockading of other inhibitory immune checkpoint receptors is being explored in hopes of providing further therapeutic options. T-cell immunoglobulin and immunoreceptor tyrosine-based inhibitory domain (TIGIT), a co-inhibitory transmembrane glycoprotein receptor expressed in various T-cell subtypes and natural killer (NK) cells, is an interesting new target for cancer immunotherapies.

“Tumor infiltrating lymphocytes (TILs) expressing TIGIT have been demonstrated in several tumor types such as non-small cell lung cancer, colorectal carcinoma, melanoma, and acute myeloid leukaemia.^{1”}


Recent studies have also shown a staining pattern correlation between TIGIT and PD-1, CD4, CD8, FOXP3 as well as other T-cell and NK cell markers.^{1,2} Biocare Medical's newly released TIGIT [BLR047F] specificity was evaluated based on multiplex IHC expressions in various T-lymphocytes, NK cells, macrophages, and dendritic cells on tonsil tissue. Co-expression level between TIGIT and utilized antibodies is shown in Table 1.

Table 1


Antibody Cocktail	Target Cells	Co-expression level*
TIGIT [BLR047F] + PD-1 [NAT105]	Activated T-cells	High
TIGIT [BLR047F] + CD3 [PS1]	T-lymphocytes	High
TIGIT [BLR047F] + CD4 [4B12]	T-helper cells	Moderate
TIGIT [BLR047F] + CD8 [C8/144B]	Cytotoxic T-cells	Low
TIGIT [BLR047F] + CD57 [NK-1]	NK cells	Moderate
TIGIT [BLR047F] + FOXP3 [86D]	Regulatory T-cells	Low
TIGIT [BLR047F] + CD68 [KP1]	Macrophages	No co-expression
TIGIT [BLR047F] + CD11c [5D11]	Dendritic cells	No co-expression

*Legend: High - >50% stained cells exhibit co-localization; Moderate – 10-49% stained cells exhibit co-localization; Low - <10% stained cells exhibit co-localization.

In this study, TIGIT [BLR047F] expressed a correct staining pattern and proved to be highly specific and sensitive to T-lymphocytes and NK cells. High level of co-localization of TIGIT and PD-1 antibodies may suggest a co-regulatory function in immune response. TIGIT plays a crucial role in inhibiting the tumor-directed immune response and is a valuable and competitive addition to the current array of antibodies for target immunotherapy research.



TIGIT (RM) [BLR047F] on tonsil



TIGIT (red) + PD-1 (brown) on tonsil

To evaluate TIGIT in your lab, contact Biocare Medical at 800-799-9499 or visit our website www.biocare.net/product/tigit-antibody/