

Probing Questions: DNA Probes vs RNA Probes vs cDNA Probes

Probing Questions: DNA Probes vs RNA Probes vs cDNA Probes

Abstract

Nucleic acid probes are invaluable tools in molecular biology, enabling clinical laboratories and researchers to detect specific DNA or RNA sequences within a sample. In this white paper, we explore the key differences between three prominent probe types: DNA probes, RNA probes, and cDNA probes. We delve into their what they are, how they are synthesized, as well as their advantages and disadvantages, guiding clinical lab staff and researchers in selecting the most appropriate probe for their specific needs.

Introduction

Nucleic acid hybridization techniques play a vital role in various clinical and research applications, including gene identification, mutation detection, viral involvement, and gene expression analysis. These techniques rely on probes – short, single-stranded fragments of DNA or RNA – that bind to complementary sequences within a target nucleic acid molecule. The specific binding allows for the detection and analysis of the target sequence.

Types of Probes

A DNA Probe is a single stranded sequence of DNA that has been synthesized from a specific area of the genome so that it can be used to search for its complimentary sequence in a sample.¹ These probes can be labeled with fluorescent tags, radioactive tags, or enzymes so that they can be visualized using microscopy. Typically, there are two ways for a DNA probe to be synthesized, using a Bacterial Artificial Chromosome (BAC) library or using DNA fragments to create a synthetic probe called an oligonucleotide.

An RNA Probe is a single stranded sequence of RNA that has been engineered to bind to a specific complimentary sequence in a sample.¹ RNA probes can target sequences of RNA OR DNA, which makes them useful in a lot of laboratory applications, like in-situ hybridization (ISH), Northern Blotting, and Southern Blotting. RNA probes typically tightly associate with its complimentary sequence in comparison to DNA probes, making them highly specific with better signal-to-noise ratios. There is only one reliable way to synthesize RNA probes through a process called in vitro transcription, which uses enzymes to copy a desired DNA sequence into a complimentary RNA nucleotide. RNA probes are also susceptible to RNase degradation, so special care is even more important and may affect their shelf life.²

A cDNA(copy DNA or complimentary DNA) probe is a single stranded DNA or RNA sequence that has been synthesized from a strand of messenger RNA (mRNA) through the use of an enzyme known as reverse transcriptase. They share the stability of DNA probes while offering the ability to target DNA or RNA sequences depending on the use case. One of cDNA's biggest advantages is the absence of introns in its sequence vs DNA and RNA probes which contain both exons and introns in their sequences. This absence of introns makes it highly specific, and typically gives a very clean final product with virtually no noise or background.

- Exons are the “coding” sections of gene sequences. They contain the instructions for building proteins and are eventually spliced together to form the final mRNA molecule that will help to create a protein. Exons can be thought of as the instructions to building a machine.
- Introns are the “non-coding” sections of gene sequences. While introns have various functions, their function isn't directly related to protein building, so they are removed during a process called RNA splicing before the mRNA is used for protein production. Consider them unnecessary information in your machine's instruction manual.

Advantages/Disadvantages

Advantages of DNA Probes:

- o **Stability:** DNA is more stable than RNA, making DNA probes easier to handle and store. They are less susceptible to degradation by enzymes like RNases, allowing for longer shelf life and less stringent handling requirements.
- o **Synthesis:** DNA probes can be readily synthesized using BAC libraries or automated DNA synthesizers, offering a faster and more cost-effective approach compared to RNA probes.
- o **Labeling:** A wider variety of labeling methods are available for DNA probes, enabling researchers to choose the most suitable detection technique for their experiment.

Disadvantages of DNA Probes:

- o **Lower Specificity:** DNA-DNA hybrids can be less stable compared to RNA-DNA hybrids. This may lead to higher background noise and reduced sensitivity, particularly when detecting targets with low abundance.
- o **Stringent Hybridization Conditions:** DNA probes often require more aggressive hybridization conditions (higher temperatures) to ensure specific binding to the target sequence. This can be detrimental to delicate tissue samples.

Advantages of RNA Probes:

- o **Higher Specificity:** RNA probes form more stable RNA-DNA hybrids compared to DNA-DNA hybrids. This translates to higher target specificity and improved signal-to-noise ratio, allowing for detection of rare transcripts or single nucleotide polymorphisms (SNPs).
- o **Lower Hybridization Temperatures:** RNA probes can often be used at lower hybridization temperatures compared to DNA probes. This is gentler on tissue samples and can help preserve their morphology.

Disadvantages of RNA Probes:

- o **Stability:** RNA is less stable than DNA, making RNA probes more susceptible to degradation by RNases. They require careful handling, storage at low temperatures, and may have a shorter shelf life.
- o **Synthesis:** Generating RNA probes typically involves in vitro transcription, which is an additional step compared to DNA probe synthesis. This can be more time-consuming, typically has a higher cost, and potentially less efficient.

Advantages of cDNA Probes:

- o **Higher Specificity:** Combine stability of DNA with targeting ability for RNA
- o **Stability:** Has the stability of DNA, don't have to worry about RNase like with RNA
- o **Copy of mRNA:** Excludes introns in its sequence, making this type of probe very clean to use

Disadvantages of cDNA Probes:

- o **Synthesis:** Requires an additional step of reverse transcription, can be time consuming to generate.

Key Differences


Feature	DNA Probes	RNA Probes	cDNA Probes
Target	DNA	DNA or RNA	DNA or RNA, dependent on use case
Synthesis	BAC Library or Oligonucleotide	In vitro transcription	Reverse transcription of mRNA
Stability	High	Moderate (RNase sensitive)	High
Specificity	Moderate	High	High
Background	Low-High	Low-Moderate	Low

Conclusion

The optimal probe selection depends on several factors, including:

- Target sequence (DNA or RNA)
- Desired level of sensitivity and specificity
- Time to synthesize or manufacture
- Stability and shelf life

DNA, RNA, and cDNA probes each offer unique advantages and disadvantages. By understanding their properties and the specific demands of the experiment, clinicians, clinical laboratories, and researchers can select the most appropriate probe to achieve optimal results in their applications.



To learn more about in situ hybridization, please visit our website at biocare.net or call 1-800-799-9499, option #3

1. J.M. Wages, NUCLEIC ACIDS I Immunoassays, Editor(s): Paul Worsfold, Alan Townshend, Colin Poole, Encyclopedia of Analytical Science (Second Edition), Elsevier, 2005, Pages 408-417, ISBN 9780123693976, <https://doi.org/10.1016/B0-12-369397-7/00730-5>

2. Robert E. Farrell, Chapter 12 - Nucleic Acid Probe Technology, Editor(s): Robert E. Farrell, RNA Methodologies (Fourth Edition), Academic Press, 2010, Pages 261-282, ISBN 9780123747273, <https://doi.org/10.1016/B978-0-12-374727-3.00012-7>

3. <https://www.genome.gov/genetics-glossary/Copy-DNA>