

Safety Data Sheet

Control Number: 907-032

Peroxidized 1 & intelliPATH™ Peroxidase Blocking Reagent
(PX968, BRR968) & (IPB5000)

June 16, 2015

BIOCARE
M E D I C A L

1. Product and Company Identification

Product(s)	Peroxidized 1 (PX968, BRR968) intelliPATH™ Peroxidase Blocking Reagent (IPB5000)
Related Product Kit	intelliPATH™ Universal HRP Detection Kit (IPK5003, IPK5011)
Intended Use & Restrictions	Immunohistochemistry reagent intended for use by professional laboratory personnel that have been properly trained in its use.
Company Information	Biocare Medical, LLC 4040 Pike Lane Concord, CA 94520 United States Tel: +1 800-799-9499 (7AM–5PM, PST, Mon-Fri) Fax: +1 925-603-8080 Web: www.biocare.net
24-Hour Emergency Contact	CHEMTREC +1 800-424-9300 (USA and Canada only) +1 703-527-3887 (all other countries)

2. Hazards Identification

Emergency Overview Clear, colorless liquid. Not particularly flammable.

Avoid contact with skin and eyes. May cause skin irritation and serious eye damage.

OSHA Hazards Target Organ Effect, Toxic by inhalation, Toxic by ingestion, Corrosive

**GHS Classification
Signal words, Symbol,
Hazard &
Precautionary
Statements**



Danger

Skin irritation (Category 2)
Serious eye damage (Category 1)

Hazard statement(s)
H315 Causes skin irritation.
H318 Causes serious eye damage.

Precautionary statement(s)
P280 Wear protective gloves/ eye protection/ face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.
Remove contact lenses, if present and easy to do. Continue rinsing.

Potential Health Effects

Skin Contact	May be harmful if exposed to skin. Causes skin burns.
Eye Contact	Causes eye burns.
Inhalation	Toxic if inhaled. Extremely destructive to the tissue of the mucous membranes and upper respiratory tract.

Safety Data Sheet

Control Number: 907-032

Peroxidized 1 & IntelliPATH™ Peroxidase Blocking Reagent
(PX968, BRR968) & (IPB5000)

June 16, 2015

BIOCARE
M E D I C A L

Ingestion Toxic if swallowed.

See Section 11: Toxicological Information for additional information.

Environmental Effects Not determined.

See Section 12: Ecological Information for additional information.

3. Composition / Information on Ingredients

Chemical Name	Synonyms	% Composition	CAS Number
Hydrogen Peroxide	N/A	<5%	7722-84-1

Note: It has been determined that the remaining ingredients of this product are not classified as hazardous according to the Federal OSHA Hazardous Communication Standard (29 CFR 1910.1200) or the Globally Harmonized System of Classification and Labeling of Chemicals.

4. First Aid Measures

General Advice Avoid further exposure to product. Consult a physician.

In case of skin contact Wash exposed area with soap and plenty of water. Remove any contaminated clothing. Seek medical attention.

In case of eye contact Wash eyes with plenty of water for at least 15 minutes. Have victim remove contact lenses. Be sure to wash under eyelids. Seek medical attention.

If inhaled Remove victim to well-ventilated area. If victim is not breathing, give artificial respiration. Seek medical attention.

If swallowed If victim is conscious, rinse mouth with water. Do not give anything to an unconscious victim. Seek medical attention.

5. Fire Fighting Measures

Suitable Extinguishing Media Water spray, alcohol-resistant foam, dry chemical, or carbon dioxide. Use extinguisher media appropriate for surrounding fire.

Special Protective Equipment for Firefighters Wear self-contained breathing apparatus and protective clothing to prevent contact with eyes and skin.

Special Exposure Hazards Avoid inhalation of fumes, which may be toxic or irritating.

6. Accidental Release Measures

Personal Precautions Use of personal protective equipment is recommended as described in Section 8. Isolate hazard area and deny entry to unprotected personnel.

Environmental Precautions Prevent further leakage or spillage if safe to do so. Keep out of drains, sewers, ditches and waterways.

Methods for containment and clean up Small spills may be contained and cleaned-up with paper towels or absorbent pads. Absorb large spills with sand or vermiculite and place in a closed container for waste disposal. Avoid physical contact during removal.

Safety Data Sheet

Control Number: 907-032

Peroxidized 1 & intelliPATH™ Peroxidase Blocking Reagent
(PX968, BRR968) & (IPB5000)

June 16, 2015

BIOCARE
M E D I C A L

7. Handling and Storage

Precautions for Safe Handling Avoid contact with eyes, skin and clothing. Do not ingest. Wash hands thoroughly after use.

Conditions for Safe Storage Keep container tightly closed. Keep container in a cool, well-ventilated area. Do not expose to excessive heat.

8. Exposure Controls / Personal Protection

Permissible Exposure Limits USA ACGIH Threshold Limit Value (Hydrogen peroxide): TWA 1 ppm
USA OSHA Table Z-1 Limits for Air Contaminants – 1910.100 (Hydrogen peroxide): 1 ppm
USA NIOSH Recommended Exposure Limits (Hydrogen peroxide): 1 ppm

Engineering Controls No special controls needed.

Personal Protective Equipment

Eye Protection Wear chemical splash goggles and face shield when eye and face contact is possible due to splashing or spraying of material.

Skin Protection Use chemically resistant gloves and a lab coat with sleeves.

Respiratory Protection Exceeding exposure limits is unlikely during normal usage. However, if irritating vapors are produced, respiratory protection is recommended. Attempt to reduce exposure levels to an acceptable range.

General Hygiene Measures Avoid contact with eyes, skin and clothing. Wash hands thoroughly after handling and before eating or drinking.

9. Physical and Chemical Properties

Appearance Clear, colorless solution

Odor None

pH (as supplied) Not determined

Melting point Not determined

Boiling point Not determined

Flash point Not determined

Flammability Not determined

Lower explosion limit Not determined

Upper explosion limit Not determined

Vapor pressure Not determined

Vapor density Not determined

Specific gravity Not determined

Solubility (in water) Not determined

Partition coefficient (n-octanol/water) Not determined

Auto-ignition temperature Not determined

Decomposition Temperature Not determined

Safety Data Sheet

Control Number: 907-032

Peroxidized 1 & IntelliPATH™ Peroxidase Blocking Reagent
(PX968, BRR968) & (IPB5000)

June 16, 2015

BIOCARE
M E D I C A L

10. Stability and Reactivity

Chemical Stability	Stable under recommended storage conditions (see Section 7).
Conditions to avoid	Keep away from incompatible materials.
Materials to avoid	Zinc, Powdered metals, Iron, Copper, Nickel, Brass, Iron and iron salts.
Hazardous decomposition products	No data available.
Possibility of hazardous reactions	No data available.

11. Toxicological Information

Acute dose effects	No data available.
Repeat dose effects	No data available.
Skin irritation	No data available.
Skin and respiratory sensitization	No data available.
Eye irritation	No data available.
Carcinogenicity	No component of this product present at levels greater than or equal to 0.1% is identified as a probable, possible, potential, known, anticipated, or confirmed carcinogen by IARC, ACGIH, NTP, or OSHA.
Mutagenicity	No data available.
Reproductive effects	No data available.
Developmental effects	No data available.
Target organ effects	Eyes, skin, respiratory system

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated at product concentrations.

12. Ecological Information

Ecotoxicity	No data available.
Persistence and degradability	No data available.
Bioaccumulation/accumulation	No data available.
Mobility in environmental media	No data available.
Other hazardous effects	No data available.

Safety Data Sheet

Control Number: 907-032

Peroxidized 1 & IntelliPATH™ Peroxidase Blocking Reagent
(PX968, BRR968) & (IPB5000)

June 16, 2015

BIOCARE
M E D I C A L

13. Disposal Considerations

Product Disposal should be in accordance with applicable national, state, and local laws and regulations. Local regulations may be more stringent than national or state requirements. Verify local and state regulations before discharging into public sewers or landfills. Do not dump into any body of water. Contact a licensed professional waste disposal service for appropriate methods of disposal.

Contaminated packaging Dispose of as unused product.

14. Transportation Information

Proper shipping name Not applicable

Technical name Not applicable

Hazard class Not applicable

Identification number Not applicable

Packing group Not applicable

DOT (US) Not regulated

IMDG Not regulated

IATA Not regulated

TDG Not regulated

15. Regulatory Information

Inventory Status

United States (TSCA) All ingredients are on the inventory or exempt from listing.

Canada (DSL / NDSL) All ingredients are on the inventory or exempt from listing.

SARA 302 Components Hydrogen peroxide (CAS 7722-84-1)

SARA 313 Components No chemicals in this product are subject to reporting requirements.

SARA 311/312 Hazards Acute health hazard, Chronic health hazard

California Prop. 65 components This product does not contain any chemicals known to the State of California to cause cancer, birth defects, or any other reproductive harm.

16. Other Information

Sources and Suitability The above information was obtained from sources available at the time of revision and believed to be accurate and reliable. The information included is not intended to be all inclusive and should only be used as a guide. Biocare shall not be held liable for any damage resulting from use, handling, or contact with the above product.

Revision Date June 16, 2015
(All SDS's before this date should be discarded.)

Replaces Revision March 8, 2013. This SDS has been updated to reflect current hazard information and product formulations.

Note: Before using this product, please visit the Biocare Medical website at <http://biocare.net> for the current version of this SDS.