

Background Terminator

Blocking Reagent

Control Number: 902-967-082317

Catalog Number: BRR 967 H, L

Description: 25, 100 ml, Ready-to-use

Intended Use:

For Research Use Only. Not for use in diagnostic procedures.

Summary & Explanation:

Background Terminator is the ultimate blocking reagent for many sources of nonspecific background staining in immunohistochemical techniques. Studies have shown it is superior to goat serum, dried milk, gelatin and BSA blocking reagents.

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Supplied As:

Goat serum plus blocking proteins in PBS with preservative and surfactant

Materials and Reagents Needed But Not Provided:

Microscope slides, positively charged
Desert Chamber* (Drying oven)
Positive and negative tissue controls
Xylene (Could be replaced with a xylene substitute*)
Ethanol or reagent alcohol
Decloaking Chamber* (Pressure cooker)
Deionized or distilled water
Wash buffer* (TBS/PBS)
Pretreatment reagents*
Enzyme digestion*
Avidin-Biotin Blocking Kit* (Labeled streptavidin kits only)
Peroxidase block*
Primary antibody*
Negative control reagents*
Detection kits*
Detection components*
Chromogens*
Hematoxylin*
Bluing reagent*
Mounting media*

* Biocare Medical Products: Refer to a Biocare Medical catalog for further information regarding catalog numbers and ordering information. Certain reagents listed above are based on specific application and detection system used.

Storage and Stability:

Store at 2°C to 8°C. Do not use after expiration date printed on vial. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user.

Staining Protocol Recommendations:

For immunohistochemistry, Background Terminator should be applied prior to primary antibody as follows:

1. Add 2-5 drops (enough to cover entire tissue section) and incubate 5-10 minutes at room temperature.
2. Rinse slides in Immunocare PBS Wash Buffer
3. Apply primary antibody.

Protocol Notes:

1. Biocare's Background Terminator is a potent blocker. In most cases, a 5 minute incubation time is adequate.
2. For severe or high background staining, incubate tissues for 15-30 minutes. Drain slides (DO NOT RINSE) and apply primary antibody.
3. Background Terminator cannot be used with two-step polymers such as MACH 3, MACH 4 or Promark detection kits.

Limitations:

This product is provided for Research Use Only (RUO) and is not for use in diagnostic procedures. Suitability for specific applications may vary and it is the responsibility of the end user to determine the appropriate application for its use.

Precautions:

1. This product is not classified as hazardous. The preservative used in this reagent is Proclin 300 and the concentration is less than 0.25%. Overexposure to Proclin 300 can cause skin and eye irritation and irritation to mucous membranes and upper respiratory tract. The concentration of Proclin 300 in this product does not meet the OSHA criteria for a hazardous substance. Wear disposable gloves when handling reagents.
2. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.
3. Microbial contamination of reagents may result in an increase in nonspecific staining.
4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.
5. Do not use reagent after the expiration date printed on the vial.
6. The SDS is available upon request and is located at <http://biocare.net>.
7. Consult OSHA, federal, state or local regulations for disposal of any toxic substances. Proclin™ is a trademark of Rohm and Haas Company, or of its subsidiaries or affiliates.

Technical Support:

Contact Biocare's Technical Support at 1-800-542-2002 for questions regarding this product.