

HPV

 | *in situ* Hybridization Probes


Results That Will Blow You Away

Biocare Medical now offers individual biotinylated probes for HPV subtypes 6, 11, 16, 18, 31 & 51


Human Papilloma Virus (HPV) is considered to be the most common sexually transmitted infection, with almost 80 million Americans infected with HPV, and about 14 million becoming newly infected each year¹. Epidemiological studies have established HPV infection as the central cause of invasive cervical cancer (ICC), with high-risk subtypes specifically associated with progression to ICC². Low-risk HPV subtypes are responsible for genital warts and a small number of recurrent respiratory papillomatosis cases^{3,4}.

- ▶ Probes available as 6 Individual ASRs – Designed to detect specific HPV subtypes 6, 11, 16, 18, 31 or 51
- ▶ Biotin labeled – Each ASR is labeled with biotin and compatible with chromogenic detection
- ▶ Choice – Choose from 6 individual ASRs for Laboratory Developed Tests


HPV Probes




HPV Type 6 Probe




HPV Type 11 Probe




HPV Type 16 Probe



HPV Type 18 Probe



HPV Type 31 Probe



HPV Type 51 Probe

Ordering Information

Product Name	Intended Use	Volume	Catalog Number
HPV Type 6 Probe	ASR	0.1 mL	BRA4030
HPV Type 11 Probe	ASR	0.1 mL	BRA4031
HPV Type 16 Probe	ASR	0.1 mL	BRA4032
HPV Type 18 Probe	ASR	0.1 mL	BRA4033
HPV Type 31 Probe	ASR	0.1 mL	BRA4034
HPV Type 51 Probe	ASR	0.1 mL	BRA4035

1. Satterwhite CL, *et al.* Sex Transm Dis. 2013 Mar;40(3):187-93 2. Clifford GM, *et al.* Br J Cancer. 2003 Jan 13;88(1):63-73 3. Garland SM, *et al.* J Infect Dis. 2009;199(6):805-14 4. Gissmann L, *et al.* Proc Natl Acad Sci USA. 1983;80(2):560-3

BIOCARE
M E D I C A L

800.799.9499
4040 Pike Lane
Concord CA 94520

www.biocare.net