

ISH


in situ Hybridization HRP Detection Kit




Spot On Results with ISH Detection

Biocare's ISH HRP Detection Kit for Biotinylated Probes is an optimized detection system incorporating Biocare's proven micro-polymer technology specifically developed to detect biotin-labeled probes on FFPE tissues. The same level of care and detail has been applied to the buffers for probe dilution and post-hybridization washes to provide the most reliable and consistent staining possible.


- ▶ Accurate – Enhanced sensitivity and specificity
- ▶ Clear – Simultaneous ISH interpretation along with assessing tissue morphology
- ▶ Distinct – Biocare's Deep Space Black™ HRP chromogen provides crisp, punctate staining
- ▶ Archivable – Chromogenic signal is stable for extended storage




Biotinylated probes targeting HPV 6 and 11



Biotinylated probe targeting HPV 16



Biotinylated probe targeting HPV 31



Biotinylated probe targeting HPV 51

Ordering Information

Product Name	Intended Use	Volume	Catalog Number
ISH HRP Detection Kit for Biotinylated Probes	IVD	6 mL	BRI4038K G
Deep Space Black™ Chromogen Kit	IVD	25 mL	BRI4015 H
DNA Hybridization Buffer	IVD	10 mL	BRI4036 G10
SSC Wash Buffer	IVD	1000 mL	BRI4039 MM
CAT Hematoxylin	IVD	500 mL	CATHE-M
Carezyme III: Pronase Kit	IVD	25 mL	PRT957K H
Pronase Buffer*	IVD	25 mL	PRB957 H

*Buffer is included in the Pronase Kit, but some users may utilize additional buffer to further dilute the enzyme for protocol optimization.

BIOCARE
M E D I C A L

800.799.9499
4040 Pike Lane
Concord CA 94520

www.biocare.net