


# SOX10<sup>(M)</sup> | *The Cleanest SOX*




## The Cleanest SOX

SOX10 (M) [BC34] from Biocare Medical is an exclusive IVD mouse monoclonal antibody that provides a highly specific and sensitive antibody for immunohistochemistry and is virtually background free.

SOX10 has been shown to stain the vast majority of cutaneous and metastatic melanomas, including 97-100% of spindle cell and desmoplastic melanomas, 100% of nevi, the majority of oligodendrogliomas and a large percentage of astrocytomas and poorly differentiated glioblastomas. PATENT PENDING.



(L) Spindle cell melanoma (DAB); (R) Pigmented melanoma (Fast red)



Melanoma stained with SOX10

Ordering Information	Concentrate	Predilute	intelliPATH™	ONCORE	VP
SOX10 (M)	ACI 3099 A, C	API 3099 AA, H	IPI 3099 G10	OAI 3099 T60	AVI 3099 G