

p40 (M)

Concentrated and Prediluted Monoclonal Antibody
902-3066-042618

BIOCARE
M E D I C A L

Catalog Number:	ACR 3066 A, C	APR 3066 AA, H	AVR 3066 KG
Description:	0.1, 1.0 ml, concentrated	6.0, 25 ml, prediluted	6.0 ml, prediluted
Dilution:	1:100	Ready-to-use	Ready-to-use
Diluent:	Van Gogh Yellow	N/A	N/A

Intended Use:

For Research Use Only. Not for use in diagnostic procedures.

Summary and Explanation:

The mouse monoclonal antibody p40 [clone BC28] recognizes an epitope unique to the p40 protein. p40 is selectively expressed in lung SCC, offering an opportunity for improved specificity (1), resulting in diminished reactivity in lung ADC and increased specificity.

The mouse monoclonal anti-p40 [BC28] demonstrated high sensitivity and specificity, staining 97% (65/67) of cases of lung SCC and 0% (0/71) of cases of lung ADC (see Performance Characteristics). p40 has also been reported in combination with TTF-1 in a method to improve specificity for SCC vs. ADC, while preserving limited tissue specimens (2,3).

Changes in expression of p40 have been implicated in other neoplastic tissues, including bladder, prostate, and head and neck cancers (1,2,3). p40 (M) [BC28] was found to be a sensitive marker in each of these tissues (see Performance Characteristics). Studies have supported the routine use of p40 as an alternative for p63 (1-4).

U.S. Patent 9,428,576 and Patents Pending.

Principle of Procedure:

Antigen detection in tissues and cells is a multi-step immunohistochemical process. The initial step binds the primary antibody to its specific epitope. A secondary antibody may be applied to bind the primary antibody, followed by an enzyme labeled polymer; or an enzyme labeled polymer may be applied directly to bind the primary antibody. The detection of the bound primary antibody is evidenced by an enzyme-mediated colorimetric reaction.

Source: Mouse monoclonal

Species Reactivity: Human; others not tested

Immunogen: a synthetic peptide corresponding to amino acids 5-17 of human p40

Clone: BC28

Isotype: IgG1

Total Protein Concentration: ~10 mg/ml. Call for lot specific Ig concentration

Epitope/Antigen: amino acids 5-17 of p40

Cellular Localization: Nuclear

Positive Tissue Control: Lung squamous cell carcinoma

Known Applications:

Immunohistochemistry (formalin-fixed paraffin-embedded tissues)

Supplied As: Buffer with protein carrier and preservative

For AVR3066KG:

p40 (M) (AVR3066G) 1 x 6ml

V-Blocker (BRR4001G) 1 x 6ml

Storage and Stability:

Store at 2°C to 8°C. Do not use after expiration date printed on vial. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

Staining Protocol Recommendations (intelliPATH and manual use):

Peroxide Block: Block for 5 minutes with Biocare's Peroxidase 1.

Pretreatment: Perform heat retrieval using Biocare's Diva Decloaker. Refer to the Diva Decloaker data sheet for specific instructions.

Protein Block (Optional): Incubate for 5-10 minutes at RT with Biocare's Background Punisher.

Primary Antibody: Incubate for 30 minutes at RT.

Probe: Incubate for 10 minutes at RT with a secondary probe.

Polymer: Incubate for 10-20 minutes at RT with a tertiary polymer.

Staining Protocol Recommendations (intelliPATH and manual use) Cont'd:

Chromogen: Incubate for 5 minutes at RT with Biocare's DAB – OR – Incubate for 5-7 minutes at RT with Biocare's Warp Red.

Counterstain:

Counterstain with hematoxylin. Rinse with deionized water. Apply Tacha's Bluing Solution for 1 minute. Rinse with deionized water.

Staining Protocol Recommendations (Ventana BenchMark XT / ULTRA Slide Staining Systems):

AVR3066 is intended for use with the Ventana BenchMark XT / ULTRA Slide Staining Systems. Refer to the User Manual for specific instructions for use. Recommended protocol parameters are as follows:

- Using **ultraView on XT / ULTRA:**

Template/Detection: ultraView DAB

Pretreatment Protocol: CC1 Mild

Primary Antibody: 32 minutes, 37°C

ultraBlock (V-Blocker BRR4001): Incubate for 4 minutes (with appropriate Option # registered by user)

V-Blocker is recommended to be applied prior to any detection system.

- Using **OptiView on ULTRA:**

Template/Detection: OptiView DAB IHC

Pretreatment Protocol: CC1 64 minutes

Peroxidase: Pre Primary Peroxidase Inhibitor

Primary Antibody: 32 minutes, 36°C

Technical Note:

This antibody, for intelliPATH and manual use, has been standardized with Biocare's MACH 4 detection system. Use TBS buffer for washing steps.

Performance Characteristics:

Nuclear staining of p40 (M) [BC28] was observed in 97% (65/67) of cases of lung squamous cell carcinoma, with no staining observed in lung adenocarcinoma cases (n=71). Staining of p40 (M) was also observed in 85.5% (41/48) of cases of urothelial carcinoma and 78% (46/59) of cases of head and neck squamous cell carcinomas. In breast cancers, only myoepithelial cells in ductal carcinoma *in situ* (DCIS) stained with p40 (M). No cases of prostate cancer were found to be positive with p40 (M). p40 (M) [BC28] nuclear staining was observed in the expected normal tissues: basal cells in prostate, myoepithelial cells in breast, urothelial cells in bladder (but not umbrella cells), stratified epithelial cells in skin, tonsil, esophagus and cervical mucosa, occasional cytotrophoblasts in placenta. (Table 2).

Limitations:

This product is provided for Research Use Only (RUO) and is not for use in diagnostic procedures. Suitability for specific applications may vary and it is the responsibility of the end user to determine the appropriate application for its use.

Precautions:

1. This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard communication and EC Directive 91/155/EC. Sodium azide (NaN₃) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for Disease Control, 1976, National Institute of Occupational Safety and Health, 1976) (5)

2. Specimens, before and after fixation, and all materials exposed to them should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens.



60 Berry Drive

Pacheco, CA 94553

USA

Rev: 062117

Tel: 800-799-9499 | www.biocare.net | Fax: 925-603-8080

p40 (M)

Concentrated and Prediluted Monoclonal Antibody
902-3066-042618

BIOCARE
M E D I C A L

Precautions Cont'd:

If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. (6)

3. Microbial contamination of reagents may result in an increase in nonspecific staining.

4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.

5. Do not use reagent after the expiration date printed on the vial.

6. The SDS is available upon request and is located at <http://biocare.net>.

Technical Support:

Contact Biocare's Technical Support at 1-800-542-2002 for questions regarding this product.

References:

1. Hibi K, *et al.* AIS is an oncogene amplified in squamous cell carcinoma. Proc Natl Acad Sci U S A. 2000 May 9; 97(10):5462-7.

2. Pelosi G, *et al.* p40 and thyroid transcription factor-1 immunoreactivity on small biopsies or cellblocks for typing non-small cell lung cancer: a novel two-hit, sparing-material approach. J Thorac Oncol. 2012 Feb; 7(2):281-90.

3. Brown AF, *et al.* Tissue-preserving antibody cocktails to differentiate primary squamous cell carcinoma, adenocarcinoma, and small cell carcinoma of lung. Arch Pathol Lab Med. 2013 Sep; 137(9):1274-81.

4. Sailer V, *et al.* Comparison of p40 and p63 expression in prostate tissues - which one is the superior diagnostic marker for basal cells? Histopathology. 2013 Jul; 63(1):50-6.

5. Center for Disease Control Manual. Guide: Safety Management, NO. CDC-22, Atlanta, GA. April 30, 1976 "Decontamination of Laboratory Sink Drains to Remove Azide Salts."

6. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition CLSI document M29-A4 Wayne, PA 2014.

VP Echelon Series antibodies are developed solely by Biocare Medical LLC and do not imply approval or endorsement of Biocare's antibodies by Ventana Medical Systems, Inc. Biocare and Ventana are not affiliated, associated or related in any way. Ventana®, BenchMark®, *ultraView* and OptiView are trademarks of Roche.

Table 1: Sensitivity and specificity of mouse monoclonal antibody p40 (M) [BC28] were determined by testing formalin-fixed, paraffin-embedded neoplastic tissues.

Pathology	Number of Specimens	Number of Positive Specimens	% Positive
Lung squamous cell carcinoma	67	65	97.0%
Lung adenocarcinoma	71	0	0%
Urothelial carcinoma	48	41	85.5%
Head and neck squamous cell carcinoma	59	46	78.0%
Breast cancer	65	18	27.6%
Prostate cancer	12	0	0%

Table 2: Tissue cross-reactivity of mouse monoclonal antibody p40 (M) [BC28] was determined by testing formalin-fixed, paraffin-embedded normal tissues.

Tissue	# positive/ total tissues	Tissue	# positive/ total tissues
Adrenal gland	0/3	Ovary	0/3
Bladder, urinary	2/3	Pancreas	0/3
Bone marrow	0/1	Parathyroid	0/3
Eye	0/1	Pituitary gland	0/2
Breast	3/3	Placenta	1/3
Brain, cerebellum	0/3	Prostate	3/3
Brain, cerebral cortex	0/3	Skin	1/1
Fallopian tube	0/3	Spinal cord	0/2
Esophagus	3/3	Spleen	0/2
Stomach	0/3	Skeletal muscle	0/3
Intestine, small intestine	0/3	Testis	0/3
Intestine, colon	0/3	Thymus	3/3
Intestine, rectum	0/3	Thyroid	0/3
Heart	0/3	Inflammatory tonsillitis*	3/3
Kidney	0/6	Ureter	3/3
Liver	0/3	Uterus cervix	3/3
Lung	0/3	Uterus (endometrium)	0/3

*B and T cells are negative. Only normal squamous epithelium is positive.