

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

1. Product and Company Identification

Product(s)	4plus Streptavidin AP Label (AP605) Background Terminator (BT967) Biotinylation Reagent (MMBR610) Da Vinci Green Diluent (PD900) Fluorescence Enhancement Probe (Mouse) (FP002) Fluorescence Enhancement Probe (Rabbit) (FP003) Fluorescence Enhancement Probe (Rabbit) (FP003) FluoroCare Anti-Fade Mountant (FP001) Horse Anti-Mouse IgG (HM606) intelliPATH Mouse Secondary Reagent (IPSC5001) intelliPATH Negative Control (ER/PR) (IPI5034) intelliPATH Universal Negative Control (IP498) MACH 1 Mouse Probe (UP537) MACH 3 Mouse Probe (MP530) MACH 3 Mouse Probe (MP537) MACH 4 Mouse Probe (UP534) MACH 4 Mouse Probe (UP536) Mouse Anti-Goat IgG (MG610) Mouse Detective (MD975) Peroxidized Diluent (PX970) Polymer Negative Control Serum (Mouse & Rabbit) (NC499) Promark Goat Probe (AP) (GP613) Promark Goat Probe (HRP) (GP626) Promark Rat Probe (RTP629) Rat and Mouse Double Stain Kit (RDS513) Renaissance Background Reducing Diluent (PD905) Universal Negative Control Serum (NC498, ORI3013)
Related Product Kit	intelliPATH Universal HRP Detection Kit (IPK5003, IPK5011) MACH 1 Universal HRP-Polymer Detection Kit (M1U539) MACH 3 Mouse AP-Polymer Detection (M3M532) MACH 3 Mouse HRP-Polymer Detection (M3M530) MACH 4 Universal AP-Polymer Kit (M4U536) MACH 4 Universal HRP-Polymer Detection (M4U534) MACH 4 Universal HRP-Polymer Detection System (BRI4012) MACH 4 Universal HRP-Polymer Kit with DAB (M4BD534) MM Biotinylation System (MMBK) Promark Goat AP-Polymer Kit (GAP514) Promark Goat HRP-Polymer Kit (GHP516) Promark Rat Detection Kit for Anti-Mouse CD31 (RT517S) Promark Rat-on-Mouse AP-Polymer Kit (RT518) Promark Rat-on-Mouse HRP-Polymer Kit (RT517)
Intended Use & Restrictions	Immunohistochemistry reagent intended for use by professional laboratory personnel that have been properly trained in its use.
Company Information	Biocare Medical, LLC 60 Berry Drive Pacheco, CA 94553 United States Tel: +1 800-799-9499 (7AM-5PM, PST, Mon-Fri) Fax: +1 925-603-8080 Web: www.biocare.net
24-Hour Emergency Contact	CHEMTREC +1 800-424-9300 (USA and Canada only) +1 703-527-3887 (all other countries)

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

2. Hazards Identification

Emergency Overview This product has been classified as non-hazardous based on the physical and/or chemical nature and/or concentration of ingredients.

Product has little to no hazards for Emergency Responders if spilled and has no unusual hazard if in a fire.

Sodium azide (<0.1%) is included as a preservative. Although it is not considered hazardous at this level, please note that accumulated sodium azide may react with lead or copper plumbing to form highly explosive metal azides. Thorough flushing of plumbing is recommended.

OSHA Hazards Not Hazardous

While this material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200), this SDS contains valuable information critical to the safe handling and proper use of the product. This SDS should be retained and available for employees and other users of this product.

GHS Classification
Signal words, Symbol,
Hazard & Not Hazardous
Precautionary
Statements

Potential Health Effects

Skin Contact May be harmful if exposed to skin. May cause skin irritation, itching, redness or inflammation.

Eye Contact May be harmful if exposed to eyes. May cause eye irritation, watering eyes, stinging or burning sensation.

Inhalation May be harmful if inhaled. May cause respiratory tract irritation, headache, dizziness, nausea or coughing.

Ingestion May be harmful if swallowed. May cause irritation of gastrointestinal tract, nausea, or vomiting.

See Section 11: Toxicological Information for additional information.

Environmental Effects None identified.
See Section 12: Ecological Information for additional information.

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

3. Composition / Information on Ingredients

Chemical Name	Synonyms	% Composition	CAS Number
Sodium azide	N/A	0-0.1%	26628-22-8
Proclin 300 (Solution of 5-Chloro-2-methyl-4-isothiazolin-3-one [2.3%] and 2-Methyl-2H - isothiazol-3-one [0.7%])	CIT, MIT	0-0.04%	26172-55-4 2682-20-4
Proclin 950 (Solution of 2-Methyl-2H - isothiazol-3-one [9.9%])	MIT	0-0.5%	2682-20-4
Proprietary Composition, multiple components	N/A	< 1% each	N/A
Water	N/A	>99%	7732-18-5

Note: It has been determined that the remaining ingredients of this product are not classified as hazardous according to the Federal OSHA Hazardous Communication Standard (29 CFR 1910.1200) or the Globally Harmonized System of Classification and Labeling of Chemicals.

4. First Aid Measures

General Advice Avoid further exposure to product. Consult a physician.

In case of skin contact Wash exposed area with soap and plenty of water. Remove any contaminated clothing. Seek medical attention.

In case of eye contact Wash eyes with plenty of water for at least 15 minutes. Have victim remove contact lenses. Be sure to wash under eyelids. Seek medical attention.

If inhaled Remove victim to well-ventilated area. If victim is not breathing, give artificial respiration. Seek medical attention.

If swallowed If victim is conscious, rinse mouth with water. Do not give anything to an unconscious victim. Seek medical attention.

5. Fire Fighting Measures

Suitable Extinguishing Media Water spray, alcohol-resistant foam, dry chemical, or carbon dioxide. Use extinguisher media appropriate for surrounding fire.

Special Protective Equipment for Firefighters Wear self-contained breathing apparatus and protective clothing to prevent contact with eyes and skin.

Special Exposure Hazards None identified.

6. Accidental Release Measures

Personal Precautions Use of personal protective equipment is recommended as described in Section 8. Isolate hazard area and deny entry to unprotected personnel.

Environmental Precautions Prevent further leakage or spillage if safe to do so. Keep out of drains, sewers, ditches and waterways.

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

Methods for containment and clean up Small spills may be contained and cleaned-up with paper towels or absorbent pads. Absorb large spills with sand or vermiculite and place in a closed container for waste disposal. Avoid physical contact during removal.

7. Handling and Storage

Precautions for Safe Handling Avoid contact with eyes, skin and clothing. Do not ingest. Wash hands thoroughly after use.

Conditions for Safe Storage Keep container tightly closed. Keep container in a cool, well-ventilated area.

8. Exposure Controls / Personal Protection

Permissible Exposure Limits	Component	Permissible Exposure Limit(s)
	Sodium Azide	ACGIH, TLV: 0.29 mg/m ³ , ceiling NIOSH, REL-C: 0.3 mg/m ³ , skin

Engineering Controls No special controls needed.

Personal Protective Equipment

Eye Protection Wear chemical splash goggles and/or face shield when eye and face contact is possible due to splashing or spraying of material.

Skin Protection Use chemically resistant gloves and a lab coat with sleeves.

Respiratory Protection Exceeding exposure limits is unlikely during normal usage. However, if irritating vapors are produced, respiratory protection is recommended. Attempt to reduce exposure levels to an acceptable range.

General Hygiene Measures Avoid contact with eyes, skin and clothing. Wash hands thoroughly after handling and before eating or drinking.

9. Physical and Chemical Properties

Appearance Clear solution. May be colorless or colored (green, red, yellow, tan), depending upon the specific product.

Odor None

pH (as supplied) Various

Melting point No data available

Boiling point No data available

Flash point No data available

Flammability No data available

Lower explosion limit No data available

Upper explosion limit No data available

Vapor pressure No data available

Vapor density No data available

Specific gravity No data available

Solubility (in water) Soluble

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

Partition coefficient (n-octanol/water)	No data available
Auto-ignition temperature	No data available
Decomposition Temperature	No data available

10. Stability and Reactivity

Chemical Stability	Stable under recommended storage conditions (see Section 7).
Conditions to avoid	Avoid buildup of sodium azide in copper or lead plumbing. Thorough flushing of plumbing with water is recommended.
Materials to avoid	Strong oxidizing agents, amines, mercaptans, reducing agents. Contact of sodium azide with heavy metals may form explosive azides. Contact of sodium azide with acids may liberate toxic gas.
Hazardous decomposition products	Hazardous decomposition products formed under fire conditions: carbon monoxide, carbon dioxide, nitrogen oxides, hydrogen chloride gas.
Possibility of hazardous reactions	Contact of sodium azide with heavy metals may form explosive azides. Contact of sodium azide with acids may liberate toxic gas.

11. Toxicological Information

Acute dose effects	No data available
Repeat dose effects	No data available
Skin irritation	Not determined
Skin and respiratory sensitization	Not determined
Eye irritation	Not determined
Carcinogenicity	No component of this product present at levels greater than or equal to 0.1% is identified as a probable, possible, potential, known, anticipated, or confirmed carcinogen by IARC, ACGIH, NTP, or OSHA.
Mutagenicity	No data available
Reproductive effects	No data available
Developmental effects	No data available
Target organ effects	No data available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated at product concentrations.

12. Ecological Information

Ecotoxicity	No data available
Persistence and degradability	No data available
Bioaccumulation/accumulation	No data available
Mobility in environmental media	No data available

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

Other hazardous effects May be harmful to the environment, particularly aquatic organisms.

13. Disposal Considerations

Product Disposal should be in accordance with applicable national, state, and local laws and regulations. Local regulations may be more stringent than national or state requirements. Verify local and state regulations before discharging into public sewers or landfills. Do not dump into any body of water. Contact a licensed professional waste disposal service for appropriate methods of disposal.

Contaminated packaging Dispose of as unused product.

14. Transportation Information

Proper shipping name Not applicable

Technical name Not applicable

Hazard class Not applicable

Identification number Not applicable

Packing group Not applicable

DOT (US) Not regulated

IMDG Not regulated

IATA Not regulated

TDG Not regulated

15. Regulatory Information

Inventory Status

United States (TSCA) All ingredients are on the inventory or exempt from listing.

Canada (DSL / NDSL) All ingredients are on the inventory or exempt from listing.

SARA 302 Components Sodium azide: CAS 26628-22-8

SARA 313 Components Sodium azide: CAS 26628-22-8, Concentration <0.1%

SARA 311/312 Hazards Sodium azide: Acute health hazard
Proclin 300: Acute health hazard
Proclin 950: Acute health hazard

California Prop. 65 components This product does not contain any chemicals known to the State of California to cause cancer, birth defects, or any other reproductive harm.

16. Other Information

Sources and Suitability The above information was obtained from sources available at the time of revision and believed to be accurate and reliable. The information included is not intended to be all inclusive and should only be used as a guide. Biocare shall not be held liable for any damage resulting from use, handling, or contact with the above product.

Revision Date January 22, 2018
(All SDS's before this date should be discarded.)

Safety Data Sheet

Control Number: 907-024

Non-Hazardous Products containing Sodium Azide & Proclin

January 22, 2018

BIOCARE
M E D I C A L

Replaces Revision April 24, 2017. This SDS has been updated to reflect current hazard information and product formulations.

Note: Before using this product, please visit the Biocare Medical website at <http://biocare.net> for the current version of this SDS.