


# intelliPATH™ Station

## Product Information

The intelliPATH Station, designed by Biocare Medical, is intended to support the intelliPATH FLX® automated slide stainer and associated components.



## Specifications

Specifications	
Weight	Approximately 275lbs
Dimensions (W x D x H)	55.1" x 28.5" x 36.3" / 140 cm x 72.4 cm x 92.3 cm


Components	
Pre-assembled workstation	
Monitor arm	

Related products (sold separately)	Cat No.
Carboy Cart	IPC100
intelliPATH FLX®	IPS0001

# Carboy Cart

## Product Information

Biocare Medical's Carboy Cart is designed to safely store and transport the IntelliPATH FLX® reagent and waste carboys. The handle enables more ergonomic movement of the cart and folds down for easy storage under the IntelliPATH™ Station or under most laboratory benches and tables.



## Specifications

Specifications	
Weight	7.2 lbs (empty) / 70 lbs loaded with full carboys
Dimensions (D x W x H)	25.5"x 13.5"x 16.5" (handle closed) / 30.5" (handle open)

Related Products (sold separately)		Cat No.
IntelliPATH™ Station	A workstation specifically designed to support the IntelliPATH FLX® instrument	IPT100