

**BIOCARE**  
M E D I C A L



## ONCORE<sup>PRO</sup>™ FISH

Dive into *in situ* hybridization (ISH) with Biocare's new automated FISH probes for ONCORE Pro. Biocare's sequence-specific probes offer superior labeling technology and are available for several cancer and disease states – yielding precise results with less ambiguity. Biocare's del-TECT™ probe design virtually eliminates error due to the truncation artifact in FFPE tissues giving you more confidence in less time.

ONCORE Pro gives you the power to automate both your IHC and FISH applications on a single instrument in a single lab. Immediately transition from running IHC to FISH using our automated PathoFISH and CytoFISH kits. Use PathoFISH on FFPE tissues and CytoFISH for urine-based samples. Run times for FISH can be completed in as little as 3.5 hours, making FISH accessible to every laboratory.



<https://biocare.net/oncore-pro>  
See Product Webpage and Datasheet for References.

**BIOCARE**  
M E D I C A L

## FISH Probes

Product Name	Probe Type	Catalog Number
PTEN del-TECT Four Color	PathoFISH	OPPR 7326 T30
TMPRSS2 /ERG del-TECT Four Color	PathoFISH	OPPR 7327 T30
ALK/EML del-TECT Four Color	PathoFISH	OPPR 7328 T30
NKX/MYC del-TECT Four Color	PathoFISH	OPPR 7329 T30
PTEN/ERG 2+2 Multiplex Four Color	PathoFISH	OPPR 7331 T30
ROS1 (6q22) Break Apart Orange/Green	PathoFISH	OPPR 7332 T30
ALK (2p23.2) Break Apart Orange/Green	PathoFISH	OPPR 7333 T30
RET (10q11.21) Break Apart Orange/Green	PathoFISH	OPPR 7334 T30
TP53 (17p13) Orange + Copy Control 17 Green	PathoFISH	OPPR 7339 T30
IGH (14q32) Break Apart Orange/Green	PathoFISH	OPRR 7340 T30
MET (7q31) Orange + Copy Control 7 Green	PathoFISH	OPPR 7341 T30
MYC (8q24) Orange + Copy Control 8 Green	PathoFISH	OPPR 7342 T30
MLL Break Apart Orange/Green	PathoFISH	OPPR 7343 T30
CytoFISH Multiplex FISH Probe	CytoFISH	OPPR 7344 T30

## FISH Kits & Ancillaries

Product Name	Probe Type	Catalog Number
ONCORE Pro FISH Kit	PathoFISH Reagents	OPRR 6064K T60
ONCORE Pro CytoFISH Kit	CytoFISH Reagents	OPRR 6068K T60
ONCORE Pro ISH Dewax Kit	FISH Reagent	OPRI 6020K T60
Wash Buffer	FISH Reagent	OPRI 6012MM
Fluoro Care Anti-Fade Mountant	FISH Reagent	FP001G5



<https://biocare.net/oncore-pro>  
See Product Webpage and Datasheet for References.

**BIOCARE**  
M E D I C A L