




BIOCARE
M E D I C A L


Key Antibodies For Dermatopathology

Top Rated Skin Antibodies




PRAME [EPR20330]

PRAME is currently being investigated as a novel diagnostic marker with ongoing clinical trials that include multiple myeloma, breast cancer, pancreatic cancer, brain cancer.




S100 Protein [15E2E2]

The S100 Protein Cocktail antibody stains melanocytes, schwannomas, peripheral neural tissue, astrocytes, benign and malignant melanomas and their metastases.



Mart-1 Cocktail

MART-1 antibody recognizes a subcellular fraction found in melanosomes. The antibody labels melanomas and tumors showing melanocytic differentiation.



Pan Melanoma + Ki-67

Pan Melanoma (MART-1 + Tyrosinase) + Ki-67 aids to identify the proliferation rate of melanocytic lesions in cases with sparse melanocytes, dense lymphocytic infiltrates or melanocytes mixed with fibroblasts.

Key Antibodies for Skin Pathology

Product Name	Source	Clone	Catalog Number
Mast Cell Tryptase [AA1]		AA1	ACI; API; OPAI; ALI; AVI; 3277
CD39 [EPR20627]		EPR20627	ACI; API; OPAI; ALI; AVI; 3283
Adipophilin		Polyclonal	ACI; API; 3138
HMB45 + MART-1 + Tyrosinase	 ⁴	HMB45 + M2-7C10/M2-9E3 + T311	CM; PM; IPI; VP 165
IFNG [IFNG/3996R]		IFNG/3996R	ACI; API 3266
Microphthalmia Transcription Factor		34CA5	CM; PM 423
Nerve Growth Factor Receptor (RM)		EP31	ACI 369
Pan Melanoma + Ki-67	 ³ 	M2-7C10/M2-9E3 + T311 + SP6	PM 362DS
pHH3 (RM)		BC37	ACI; API 3130
SOX10		BC34	ACI; API; IPI; AVI 3099
CD34		QBEnd/10	CM; PM; IP 084
Cytokeratin [AE1] LMW		AE1	PM 081
Factor XIIIa		E980.1	CM; PM; IP 357
HMB45		HMB45	CM; PM; IP 057
MART-1 + Tyrosinase + pHH3	 	BC37, M2-7C10 + M2-9E3, T311	API 3168
MART-1 + Tyrosinase + SOX10		BC34,M2-7C10 + M2-9E3, T311	API 3165
Mart-1 Cocktail		M2-7C10 + M2-9E3,	CM; PM; IP 077
Melanoma Cocktail		HMB45 + M2-7C10 + M2-9E3	CM; PM; VP 078
p16 + PRAME	 	EPR20330 + BC42	API 3256DS
PRAME [EPR20330]		EPR20330	ACI; API; OPAI; ALI; AVI 3252
S100 Protein [15E2E2] (M)		15E2E2	CM; PM 128
Melan A		A103	ACI; API; 3114
MLH-1		G168-15	CM; PM; IPI; AVI 220
MSH2		FE11	CM; PM; AVI 219
MSH6		44	CM; PM; IPI; AVI 265
PMS2		A16-4	CM; PM; IPI; AVI 344
MLH-1 [BC23]		BC23	ACI; API 3214
MSH6 [BC19]		BC19	ACI; API 3215
HHV-8		13B10	ACA; APA; OPAA; AVA; ALA 3290